

Sterile Storage Carts

Features:

- Hook configuration to accommodate 3 – 4 baskets
- Hooks tested to 50kg load capacity
- 4x 75mm non marking swivel castors (2x braked).

Dimensions:

Cart:

Width (mm) 610, Depth (mm) 470, Height (mm) 1450

Individual Basket:

Width (mm) 575, Depth (mm) 280, Height (mm) 265

Ordering:

Product	Code
Single Sterile Storage Cart	STC300
Double Sterile Storage Cart	STC600

Hygienic solutions for your storage.

Our single and double-sided Storage Carts are made from 304 grade strong stainless steel with a seamless welded frame for durability and easy cleaning. Lockable swivel castors provide stability, while the hook configuration offers versatile options to accommodate baskets.

Options:

Single Bowl Stand



Double Bowl Stand

