

Surgical Steps

A solid step up.

Our stable 4-footed steps are seamlessly welded from 304 grade stainless steel with removable stainless steel ribbed tread for easy cleaning. Non-slip, non-marking, adjustable rubber feet provide additional stability and reliability.

Options:



Features:

- 2 non-slip, non-marking adjustable rubber feet
- 2 conductive polymer screw-in feet
- Load capacity: 150kg.

Dimensions:

Single Step (130mm High):

Width (mm) 600, Depth (mm) 350, Height (mm) 130

Single Step (220mm High):

Width (mm) 600, Depth (mm) 350, Height (mm) 220

Double Step:

Width (mm) 600, Depth (mm) 650, Height (mm) 220/440

Ordering:

Product	Code
Single Step (130mm High)	STE1631
Single Step (220mm High)	STE1632
Double Step	STE26624