

Theatre Storage Trolleys

Specifically designed for hospital theatres.

The Agile Medical MTC Enclosed Transport Trolley is a specially designed, high capacity mobile storage cabinet with high impact resistance, easy clean qualities and a Stainless Steel body that is impervious to chemicals and high energy radiation such as X-rays and Gamma rays.

Features:

- 8x 100mm ABS Trays, 1x 200mm ABS Tray
- Easily manoeuvrable
- Lockable castors.

Dimensions:

Width (mm) 740, Depth (mm) 580, Height (mm) 1655

Ordering:

Product	Code
Theatre Storage Trolley - Drawers	MTC1622
Theatre Storage Trolley - Solid Doors	MTC1676

