

# Three Section Couches

For a wider range of uses.

Multi-discipline three-section models suit a variety of medical, clinical and minor surgical applications. Electronic variable height is included as standard, with selected couches available with a hydraulic lift alternative.



Retractable wheel base with foot lever control allows couch to be moved/positioned with ease.

## Attachments to order:



# Three Section Couches

## Features:

- Gas assisted adjustable backrest and footrest sections
- Retractable wheelbase with foot lever control
- Anti-microbial seamless 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating.

## Dimensions:

**Length:** 1900mm **Height adjustment:** 480-890mm **Safe working load:** 225kg

## Vinyl Colours:

<b>1</b> Burnt Red	<b>2</b> Rose	<b>3</b> Mulberry	<b>4</b> Lilac	<b>5</b> Navy
<b>6</b> Royal	<b>7</b> Lido	<b>8</b> Aqua	<b>9</b> Ocean	<b>10</b> Citrus
<b>11</b> Sage	<b>12</b> Grey	<b>13</b> Dove	<b>14</b> Ginger	<b>15</b> Buttercup
<b>16</b> Taupe	<b>17</b> Barley	<b>18</b> Stone	<b>19</b> White	<b>20</b> Black

## Options to order:



## Ordering:

Add the chosen colour code to the end of the product code.

For example: **CC3MD/25H6** for Royal Three Section Couch - Hydraulic

Product	Code	
Three Section Couch - Electronic	<b>CC3MD/25E</b>	
Three Section Couch - Hydraulic	<b>CC3MD/25H</b>	
Attachments to order		
Hygiene roll holder	<b>HRH001</b>	
Safety side rails (pair)	<b>SSR001</b>	
Safety side pads (pair)	<b>SSP001</b>	
Head cushion support	<b>HCS001</b>	
IV pole	<b>IVP001</b>	
Ankle stirrups (pair)	<b>ANS001</b>	
Leg troughs (pair)	<b>LGT001</b>	
Phlebotomy arm (single / pair)	<b>PAS001</b>	<b>PAP001</b>
Adjustable arm supports (pair)	<b>AAS001</b>	
Options to order		
Additional +3" width	<b>A3WB</b>	
Breather hole with vinyl wedge to plug hole	<b>BHW001</b>	
Additional top layer of anatomic memory foam	<b>AMF001</b>	
Alternative electronic footswitch	<b>AEF001</b>	
Base variant to allow hoist access	<b>BVH001</b>	