

Two Section Couches

Multi-discipline two-section models suit a variety of medical, clinical and minor surgical applications. Electronic variable height comes as standard, with an optional hydraulic lift alternative available on selected couches.



Retractable wheel base with foot lever control allows couch to be moved/positioned with ease.

Attachments to order:



*CC2MD/25 & CC2MD/28 only

Two Section Couches

Features:

- Gas assisted adjustable backrest section
- Retractable wheelbase with foot lever control
- Anti-microbial seamless 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating.

Dimensions:

Length: 1900mm **Height adjustment:** 480-890mm **Safe working load:** 225kg

Vinyl Colours:

1 Burnt Red	2 Rose	3 Mulberry	4 Lilac	5 Navy
6 Royal	7 Lido	8 Aqua	9 Ocean	10 Citrus
11 Sage	12 Grey	13 Dove	14 Ginger	15 Buttercup
16 Taupe	17 Barley	18 Stone	19 White	20 Black

Options to order:



Ordering:

Add the chosen colour code to the end of the product code.

For example: **CC2MD/25H5** for Rose Two Section Couch - Hydraulic

Product	Code	
Two Section Couch - Electronic	CC2MD/25E	
Two Section Couch - Hydraulic	CC2MD/25H	
Attachments to order		
Hygiene roll holder	HRH001	
Safety side rails (pair)	SSR001	
Safety side pads (pair)	SSP001	
Head cushion support	HCS001	
IV pole	IVP001	
Ankle stirrups (pair)	ANS001	
Leg troughs (pair)	LGT001	
Phlebotomy arm (single / pair)	PAS001	PAP001
Adjustable arm supports (pair)	AAS001	
Options to order		
Additional +3" width	A3WB	
Breather hole with vinyl wedge to plug hole	BHW001	
Additional top layer of anatomic memory foam	AMF001	
Alternative electronic footswitch	AEF001	
Base variant to allow hoist access	BVH001	