

# Ultrasound Scanning/Medical Procedure Chairs

## Precise patient positioning.

Choose from electronic height operation with foot console or hydraulic height operation with pneumatic foot pump to achieve optimal patient positioning. The long backrest is ideal for ultrasound examinations, while the compact base allows close access to the patient. The steel frame and vinyl upholstery both benefit from an anti-microbial protective coating for maximum hygiene.



Provides precise patient positioning with raised seat tilt facility.

### Attachments to order:



HCS001



LGT001



AAS001

### Options to order:



AMF001

# Ultrasound Scanning/Medical Procedure Chairs

## Features:

- Variable height
- Gas assisted adjustable back and mid sections
- Manually adjusted leg section
- Seat section has hydraulic tilt of 20°
- Retractable wheelbase with foot lever control.

## Dimensions:

**Length:** 1860mm **Height adjustment:** 580-860mm **Safe working load:** 225kg

## Vinyl Colours:



## Ordering:

Add the chosen colour code to the end of the product code.

For example: **CH3PR/24E16** for Taupe Ultrasound / Scanning / Clinic Chair - Electronic

Product	Code
Ultrasound / Scanning / Clinic Chair - Electronic	<b>CH3PR/24E</b>
Ultrasound / Scanning / Clinic Chair - Hydraulic	<b>CH3PR/24H</b>
Attachments to order	
Head cushion support	<b>HCS001</b>
Leg support troughs (pair)	<b>LGT001</b>
Adjustable arm supports (pair)	<b>AAS001</b>
Options to order	
Additional top layer of anatomic memory foam	<b>AMF001</b>