

Under Worktop Procedure Trolleys

Low height for high performance.

Designed to fit under standard worktops, the Low Height Procedure Trolley enables maximum utilisation of space for storing essential equipment. Integrated push handle and an extendable work surface increases versatility versus fixed units, while lockable drawers provide increased security and hygiene.



Options:

Narrow Trolley - Under Worktop



520mm

LX Trolley - Under Worktop



690mm

Double Trolley - Under Worktop



1125mm

Under Worktop Procedure Trolleys

Features:

- Durable epoxy-coated metal construction with 10-year warranty
- One piece, fluid retaining ABS top for easy cleaning.
- Protective corner profiles and base bumper.

Colours:



Dimensions:

Narrow Trolley:

Width (mm) 520, Depth (mm) 540, Height (mm) 850

LX Trolley:

Width (mm) 690, Depth (mm) 540, Height (mm) 850

Double Trolley:

Width (mm) 1125, Depth (mm) 545, Height (mm) 850

Ordering:

Product	Drawers	Col	Code
Narrow Procedure Trolley - Under Worktop			NC28041
			NC28042
LX Procedure Trolley - Under Worktop			LX28PRO1
			LX28PRO2
Double Procedure Trolley - Under Worktop			DC34091
			DC34092
Accessories			
Magnetic label holders			SUPPLIED
Drawer divider set (2 supplied, additional available)			SUPPLIED
IV pole			IVP100
Sharps bin holder			SBH100
Catheter holder			CHS114
1 box glove holder			GDS119
Gel holder			GMH100
Security clamp			DRS200
Waste container with lid			WCS213
Multi-storage basket			MSL100
File holder			MSL100

For our full range of trolley accessories see pages **108-111**.