

# Vena Vario Phlebotomy Chairs

## Combine comfort and seamless functionality.

Easy operation makes it possible to act quickly in emergency situations, with the patient able to be moved safely and quickly into the Trendelenburg position. The Vario Levant has 2 x 'Vena® fixed' with PU arm pads, integrated hydraulic foot pump operation and features a synchronously adjustable back and leg rest allowing the operator to place the patient in the optimum position when carrying out procedures.

## Options:



## Features:

- Available in Stamskin® and PUxx Material options
- Wipe clean hygienic design.
- Robust steel base, white (RAL 9016) with 4 feet
- Protective glides
- Medical device in accordance with MDR 2017/745 (class I).

## Vario Levant:

- Hydraulic height adjustment 540-740mm
- Hand controlled synchronic back and leg rest operation.

## Stamskin® Colours:

- |   |   |   |   |
|---|---|---|---|
|  K85 Anthracite  |  K05 Black         |  K294 Blush        |  K170 Burgundy Red |
|  K189 Capri Blue |  K89 Dark Blue     |  K196 Dutch Orange |  K292 Forest Green |
|  K295 Mouse Grey |  K190 Ochre Yellow |  K394 Pink Orchid  |  K258 Taupe        |

## Ordering:

Add the chosen colour code to the end of the product code.  
For example: **PC2001/MAK29** for Dark Blue Vena Vario Phlebotomy Chair.

| Product                            | Code             |
|------------------------------------|------------------|
| Vena Vario Phlebotomy Chair        | <b>PC2001/MA</b> |
| Vena Vario Levant Phlebotomy Chair | <b>PCH2001</b>   |