

infinity

CSSD

Static - Single Bay

Product	Weight	Code
 <p>Static, Single Bay - 2 x Wire Grid Shelves, 5 x Telescopic Shelves</p>	30.2kg	MDXST300-1-F/2W5T
 <p>Static, Single Bay - 8 x Wire Grid Shelves</p>	21.4kg	MDXST300-1-F/8W

Mobile - Single Bay

Product	Weight	Code
 <p>Mobile, Single Bay - 2 x Wire Grid Shelves, 5 x Telescopic Shelves</p>	31.1kg	MDXST300-1-W/2W5T
 <p>Mobile, Single Bay - 8 x Wire Grid Shelves</p>	22.3kg	MDXST300-1-W/8W

Static - Double Bay

Product	Weight	Code
 <p>Static, Double Bay, 4 x Wire Grid Shelves & 10 x Telescopic Shelves</p>	56.9kg	MDXST300-2-F/4W10T
 <p>Static, Double Bay, 16 x Wire Grid Shelves</p>	38.0kg	MDXST300-2-F/16W

Mobile - Double Bay

Product	Weight	Code
 <p>Mobile, Double Bay, 4 x Wire Grid Shelves & 10 x Telescopic Shelves</p>	58.2kg	MDXST300-2-W/4W10T
 <p>Mobile, Double Bay, 16 x Wire Grid Shelves</p>	39.3kg	MDXST300-2-W/16W

Designed for flexibility.

Agile Medical's innovative, high-quality aluminium product is a flexible system that expands easily to different widths and depths.

It is specifically designed for safe and protected storage of sterile containers, with options of wire shelves for lightweight items and pull-out support runners for heavier items.

Key benefits:

- Protective anodized aluminium construction
- Adjustable pull-out shelves to accommodate all sizes of sterile container.

Dimensions:

Static Single Bay:

Width (mm) 394, Depth (mm) 642, Height (mm) 1940

Static Double Bay:

Width (mm) 737, Depth (mm) 642, Height (mm) 1940

Mobile Single Bay:

Width (mm) 394, Depth (mm) 642, Height (mm) 1996

Mobile Double Bay:

Width (mm) 737, Depth (mm) 642, Height (mm) 1996

